

# EDGE 3 *Stretto*



The Q6 Edge 3 Stretto features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12,5" or 14" drive wheels and optional 305 mm of iLevel power-adjustable seat height.



**QUANTUM®**

#1 FOR REHAB POWER



| [quantumrehab.eu](http://quantumrehab.eu)

## Specifications

DRIVE WHEELS	14" or 12,5" pneumatic or flat free
CASTOR WHEELS	Front: 5" flat free Rear: 6" flat free
MAXIMUM SPEED <sup>1</sup>	Up to 10 km/h Up to 7,2 km/h w/iLevel activated
GROUND CLEARANCE	74 mm
TURNING RADIUS <sup>2</sup>	533 mm
OVERALL LENGTH <sup>2</sup>	826 mm
BASE WIDTH	530 mm with 12,5" drive wheels 552 mm with 14" drive wheels
<b>SEATING SIZES</b>	
Synergy TRU-Balance 3 Static Seat	Width: 305 mm - 508 mm Depth: 305 mm - 508 mm
TRU-Balance 3 Power Tilt	Width: 305 mm - 508 mm Depth: 305 mm - 508 mm
TRU-Balance 3 Tilt and iLevel	Width: 305 mm - 508 mm Depth: 305 mm - 508 mm
TRU-Balance Flex Seating	Width: 381 mm - 508 mm Depth: 381 mm - 508 mm
<b>SEAT-TO-FLOOR HEIGHT</b>	
Synergy TRU-Balance 3 Static Seat <sup>5</sup>	406 mm - 457 mm
TRU-Balance 3 Power Tilt <sup>5</sup>	445 mm - 495 mm
TRU-Balance 3 Tilt & iLevel <sup>5</sup>	445 mm - 495 mm
TRU-Balance Flex Seating	419 mm - 464 mm
<b>MANUAL RECLINE <sup>6</sup></b>	
	Yes (TRU-Balance Flex only)
<b>BATTERY SIZE <sup>3</sup></b>	
	2x 40 or 50 Ah
<b>BATTERY WEIGHT <sup>7</sup></b>	
	15 kg each
<b>ELECTRONICS <sup>8</sup></b>	
	Q-Logic 3 or NE+
<b>BATTERY CHARGER</b>	
	8 Amp, off-board
<b>MOTOR PACKAGES</b>	
	2-Pole, Mid-Wheel <sup>6</sup>
<b>WEIGHT CAPACITY <sup>4</sup></b>	
	136 kg
<b>BASE WEIGHT</b>	
	57,2 kg (without batteries)

- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- Without front riggings
- Battery capacity may change when adding power positioning systems.
- Weight capacity is dictated by the selection of seating.
- All measurements +/- 3 mm with 457 x 457 mm seat, fixed towers, 136 kg weight capacity and solid tires (add 10 mm with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- Battery weight may vary +/- 1,3 kg based on manufacturer.
- All Quantum Rehab units can be ordered with specialty controls.

## Features

- Independent SRS (Smooth Ride Suspension)
- iLevel option available with up to 305 mm of seat elevation at walking speed (up to 7,2 km/h)
- Fender lights and USB charger included



## Colours\*

\* See website for a complete list of colour options.



Cherry Bomb

Blue Jean

Lemon Crush

Raspberry Beret

Evolution Orange



Sugar Plum

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



**Pride Mobility Products BV**  
De Zwaan 3  
1601 MS Enkhuizen  
The Netherlands

Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039  
quantumrehab.eu · info@prideeurope.nl



EU\_SSQ0005/REVA/01DEC2020