

EDGE 3 *Stretto*



The **Q6 Edge 3 Stretto** features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12,5" or 14" drive wheels and optional 305 mm of iLevel power-adjustable seat height.



Extreme narrow width
of just 520 mm with
12,5" drive wheels



QUANTUM[®]

#1 FOR REHAB POWER



| quantumrehab.eu

Specifications

DRIVE WHEELS	14" or 12,5" pneumatic or flat free
CASTOR WHEELS	Front: 5" flat free Rear: 6" flat free
MAXIMUM SPEED ¹	Up to 10 km/h
GROUND CLEARANCE	74 mm
TURNING RADIUS ²	533 mm
OVERALL LENGTH ²	826 mm
BASE WIDTH	530 mm with 12,5" drive wheels 552 mm with 14" drive wheels
SEATING SIZES	
TRU-BALANCE FLEX	Width: 406 - 508 mm Depth: 406 - 508 mm
TRU-BALANCE 3	Width: 305 - 610 mm Depth: 305 - 610 mm
SEAT TO FLOOR HEIGHT ^{5,6,7}	420 mm (Low) 445 mm (Middle) 470 mm (High)
MANUAL RECLINE	Yes (TRU-Balance Flex only)
BATTERY SIZE ³	Two 12 volt, deep-cycle 40 or 62 Ah, AGM
BATTERY WEIGHT ^{5,6}	40 Ah: 12,5 kg each 62 Ah: 18,7 kg each
ELECTRONICS	Q-Logic 3, Q-Logic 3e or NE+
BATTERY CHARGER	8 Amp, off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6
WEIGHT CAPACITY ³	136 kg
BASE WEIGHT	57,2 kg (without batteries)

1. Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure
2. Without front riggings
3. Capacity may change when adding power positioning systems
4. Battery weight may vary +/- 1,4 kg based on manufacturer
5. Seat-to-floor heights are measured with the unit unloaded
6. Tilt or lift adds 35 mm
7. Measurements without seat cushion

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel option available with up to 305 mm of seat elevation at walking speed (up to 7,2 km/h)
- Extreme narrow width of 520 mm with 12,5" drive wheels, 552 mm with 14" drive wheels
- Midwheel drive
- Speed 10 km/h
- 12,5" pneumatic drive wheels
- Max user weight capacity 136 kg
- Fender lights and USB charger optional



Colours



Pride Mobility Products BV
De Zwaan 3
1601 MS Enkhuizen
The Netherlands

Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039
quantumrehab.eu · info@prideeurope.nl



EU_SSQ0005/REVB/01APR22