

EDGE 3



The **Edge 3** has upgraded SRS (Smooth Ride Suspension) providing a smooth, comfortable ride and advanced stability. Optional iLevel is available at 7,2 km/h with up to 305 mm of power adjustable seat height. Standard features include LED fender lights, a USB charger, and drive wheel colour accents.



Now available with
optional 11,3 km/h
4-pole motor



QUANTUM®

#1 FOR REHAB POWER



| quantumrehab.eu

Specifications

DRIVE WHEELS	14" solid or pneumatic
CASTER WHEELS	6" solid
MAXIMUM SPEED ¹	Standard: Up to 10 km/h ¹⁰ Optional: Up to 11,3 km/h ¹¹ Up to 7,2 km/h with iLevel
GROUND CLEARANCE	69 mm
TURNING RADIUS ²	521 mm
LENGTH ²	820 mm
WIDTH	622 mm
SEATING SIZES	
Synergy TRU-Balance 3 Static Seat	Width: 254 - 559 mm Depth: 254 - 559 mm
TRU-Balance 3 Power Tilt	Width: 305 - 559 mm Depth: 305 - 559 mm
TRU-Balance 3 Tilt and iLevel	Width: 305 - 559 mm Depth: 305 - 559 mm
TRU-Balance Flex Seating	Width: 381 - 559 mm Depth: 305 - 508 mm
SEAT TO FLOOR HEIGHT ⁹	
Synergy/TRU-Balance 3 Static Seat ⁷	154 - 440 mm
TRU-Balance 3 Power Tilt ⁷	445 - 495 mm
TRU-Balance 3 Lift & Tilt ⁷	451 - 512 mm
TRU-Balance Flex Seating	419 - 495 mm
MANUAL RECLINE ⁵	Yes
BATTERY SIZE ³	2x NF-22 (50, 55 or 62 Ah)
BATTERY WEIGHT ⁸	18 kg each
AVAILABLE ELECTRONICS ⁶	Q-Logic 3
BATTERY CHARGER	8 Amp, off-board
DRIVETRAIN	2-Pole, Mid-Wheel 6® 4-Pole, Mid-Wheel 6®
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT	63 kg

- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- Without front riggings.
- Battery capacity may change when adding power positioning systems.
- Weight capacity is dictated by the selection of seating.
- Maximum dimensions of the Manual Recline vary by power base and seating system. Contact us for details.
- All Quantum Rehab units can be ordered with specialty controls. Contact us for availability.
- All measurements +/- 3 mm with 460 x 460 mm seat, fixed towers, 136 kg weight capacity and solid tires (add 9 mm with pneumatic tires). Widths may increase seat-to-floor height. Contact customer service for details.
- Battery weight may vary +/- 1,4 kg based on manufacturer.
- 419 mm seat-to-floor height is not available with power tilt or iLevel. 495 mm seat-to-floor height is not available with static seat or power recline.
- 10 km/h speed available only with Q-Logic 3.
- Up to 11,3 km/h speed available only with 4-pole motor option and Q-Logic 3 or Q-Logic 3e.

NOTE: All specifications subject to change without notice.


Features

- Upgraded SRS (Smooth Ride Suspension)
- 7,2 km/h with 305 mm iLevel
- Fender lights and USB charger included
- Drive wheel colour accents included



Colours*

* See website for a complete list of colour options.



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



Pride Mobility Products BV
De Zwaan 3
1601 MS Enkhuizen
The Netherlands

Tel: +31 (0)85-0408630 • Fax: +31 (0)85-0430039
quantumrehab.eu • info@prideeurope.nl



EU_SSQ0004/REVB/09NOV2021