

# Q6 edge<sup>®</sup> Z

---





The **Q6 Edge Z** utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6 design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability.



**QUANTUM<sup>®</sup>**

#1 FOR REHAB POWER

  | [quantumrehab.eu](http://quantumrehab.eu)

# Specifications

<b>DRIVE WHEELS</b>	14" pneumatic or solid
<b>CASTER WHEELS</b>	6" solid
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 10 km/h
<b>GROUND CLEARANCE <sup>2</sup></b>	76 mm
<b>TURNING RADIUS <sup>2</sup></b>	502 mm
<b>LENGTH <sup>2,4</sup></b>	904 mm
<b>WIDTH <sup>2,4</sup></b>	667 mm
<b>SEATING SIZES</b>	
<b>TRU-Balance 3 Static</b>	Width: 305 - 610 mm Depth: 305 - 610 mm
<b>TRU-Balance 3 Power Lift</b>	Width: 305 - 610 mm Depth: 356 - 610 mm
<b>TRU-Balance 3 Lift</b>	Width: 305 - 610 mm Depth: 305 - 610 mm
<b>SEAT-TO-FLOOR HEIGHT <sup>5</sup></b>	
<b>TRU-Balance 3 Static</b>	419 - 470 mm
<b>TRU-Balance 3 Power Lift</b>	470 - 521 mm
<b>TRU-Balance 3 Lift</b>	470 - 521 mm
<b>BATTERIES <sup>3</sup></b>	
<b>BATTERY CHARGER</b>	2x 55, 62 or 75 Ah
<b>BATTERY CHARGER</b>	8 Amp, off-board
<b>ELECTRONICS</b>	Q-Logic
<b>MOTOR PACKAGES</b>	4-Pole Mid-wheel-6 drive
<b>WEIGHT CAPACITY</b>	136 kg
<b>BASE WEIGHT</b>	65 kg

1. Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure
2. Without front riggings
3. Battery capacity may change when adding power positioning systems
4. All measurements +/- 3 mm with 457 x 457 mm seat, fixed towers, 136 kg weight capacity and solid tires (add 10 mm with pneumatic tires). Depths may increase seat-to-floor height
5. Listed measurements are in the lowest position; does not include elevating system unless specified

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for base model only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

# Features

- Mid-Wheel 6 drive allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance 3 Power Positioning Systems
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- Standard switch-operated LED fender lights

# Options

- Power adjustable seat height
- TRU-Balance 3 Seating
- Q-Logic Drive Controls
- 12,8 km/h motor package
- Accu-Trac Advanced
- Tracking Technology
- Trilogy vent tray
- Compact vent tray

# Colours



**Pride Mobility Products BV**  
De Zwaan 3  
1601 MS Enkhuizen  
The Netherlands  
Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039  
quantumrehab.eu · info@prideeurope.nl



EU\_SSQ0007/REVC/020FEB2022