

Q6 edge[®] Z



The **Q6 Edge Z** utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6 design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability.



QUANTUM[®]

#1 FOR REHAB POWER



| quantumrehab.eu

Specifications

DRIVE WHEELS	14" pneumatic or flat free, also available as extra wide high profile tires
CASTER WHEELS	9" flat free
MAXIMUM SPEED ¹	Up to 10 km/h 12,8 km/h optional
GROUND CLEARANCE ²	76 mm
TURNING RADIUS ²	502 mm
LENGTH ^{2,4}	904 mm
WIDTH ^{2,4}	667 mm
SEATING SIZES	
SEATING SIZES ⁵	Width: 305 - 610 mm Depth: 305 - 610 mm
SEAT TO FLOOR HEIGHT ^{5,8,9,10}	420 mm (Low) 445 mm (Middle) 470 mm (High)
BATTERIES ³	Two 12 volt, deep-cycle 62 or 75 Ah
BATTERY WEIGHT	62 Ah: 18,7 kg 75 Ah: 24,4 kg
BATTERY CHARGER	8 Amp, off-board
ELECTRONICS	Q-Logic
MOTOR PACKAGES	4-Pole Mid-wheel-6 drive
WEIGHT CAPACITY	136 kg
BASE WEIGHT	65 kg

- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure
- Without front riggings
- Battery capacity may change when adding power positioning systems
- All measurements +/- 3 mm with 457 x 457 mm seat, fixed towers, 136 kg weight capacity and solid tires (add 10 mm with pneumatic tires). Depths may increase seat-to-floor height
- Listed measurements are in the lowest position; does not include elevating system unless specified

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for base model only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Features

- Mid-Wheel 6 drive allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance 3 Power Positioning Systems
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- Standard switch-operated LED fender lights

Kleuren



Red
Corvette



Electric
Blue



Black
in
Black



Green
Machine



Night
Train



Grey
Street



Deep
Purple



Evolution
Orange



Pretty in
Pink



Midnight
Blue



Pride Mobility Products BV
De Zwaan 3
1601 MS Enkhuizen
The Netherlands

Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039
quantumrehab.eu · info@prideeurope.nl



EU_SSQ0007/REVD/01APR22