



MOTABILITY ID 706910

EDGE 3 STRETTO - CONTOURED SEAT / 6 MPH

DEALERS DETAILS

Dealers Name _____
 Order Date _____
 Address _____
 _____ Postcode _____
 Phone _____ Email _____
 PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____
 _____ Postcode _____
 Quantum Specialist _____
 ⓘ Clients Weight (required for suspension settings) _____

BASE MODEL

- 14 Colour Options
- Flip-up Armrests
- 70° Swing-Away Legrests
- Contoured Seat with Fold-Forward Backrest
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 6 Mph (9.5 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	18"W x 18"D	Standard	20"W x 20"D	FOC
-------------	-----	-------------	----------	-------------	-----

Seat package includes standard headrest and auto buckle lap belt.

SEAT CUSHIONS

Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Seat cushion dimensions will match seat package selected above.

SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard		
Manual Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard			
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>	FOC			
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	FOC			
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range 10" to 16". 5"W x 11"L Footplates.</i>	FOC
Seat pan to heel measurement	(Leave blank for standard factory)

NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>	Standard
VR2 Controls. <i>With no power options.</i>	FOC

HAND CONTROL MOUNTS

Fixed.	Standard			
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>	FOC			
Mounting Position.	Left	FOC	Right	Standard

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.

Standard

NOTES & SPECIAL REQUESTS

3



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
Motability ID 706910 - 17/03/23 V1



Crash Tested



V1: 14/06/2022

