



MOTABILITY ID 721039

Q4 WITH TB FLEX SEAT / POWERED TILT / POWER RECLINE

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

- 8 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors



- ATX Suspension System
- Powered Tilt
- Powered Recline
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	White Satin	FOC
Midnight Blue	FOC	Red Corvette	FOC
Electric Blue	FOC	Please download our free PQ-365 app to view colours and other product features in more detail.	
Deep Purple	FOC		
Pretty in Pink	FOC		
Evolution Orange	FOC		



SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	Standard
Seat to Floor Height. <i>From seat pan.</i>	16" 17" <i>Changes Tilt range to -5° to 25°</i> 18" <i>Changes Tilt range to -10° to 20°</i>

STANDARD BACK SUPPORT POSITIONING

Powered Back Recline. <i>Adjustment range 95° to 125°. Uses 1 power option.</i>	Standard
---	----------

STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
---	----------

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	FOC
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	FOC
Arm Support Pads	10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range 10" to 16". 5"W x 11"L Footplates.</i>	FOC

Seat pan to heel measurement (Leave blank for standard factory) Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard (STD)
--	----------------

NON-EXPANDABLE CONTROLS

NE + Controls. <i>2 power option ready.</i>	Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>
	Attendant Control. <i>Hand control mounted at rear.</i>
	FOC

HAND CONTROL MOUNTS

Fixed.	Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>	FOC
Mounting Position.	Left FOC Right Standard

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 721039 - 18/01/23 V1



Crash Tested



V1: 14/06/2022

