



MOTABILITY ID 1000791

## Q4 WITH TB FLEX SEAT / PWR TILT / PWR RECLINE / PRO BACK

### DEALERS DETAILS

Dealers Name \_\_\_\_\_  
 Order Date \_\_\_\_\_  
 Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Phone \_\_\_\_\_ Email \_\_\_\_\_  
 PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Quantum Specialist \_\_\_\_\_  
 Clients Weight \_\_\_\_\_

### BASE MODEL

- 7 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors



- ATX Suspension System
- Powered Tilt
- Powered Recline
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

### COLOURS

Back in Black	Standard	Evolution Orange	FOC
Midnight Blue	FOC	Red Corvette	FOC
Electric Blue	FOC	Please download our free PQ-365 app to view colours and other product features in more detail.	
Deep Purple	FOC		
Pretty in Pink	FOC		



### SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

### SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

## SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>				Standard
Seat to Floor Height. <i>From seat pan.</i>	16"	17" <i>Changes Tilt range to -5° to 25°</i>	18" <i>Changes Tilt range to -10° to 20°</i>	

## STANDARD BACK SUPPORT POSITIONING

Powered Back Recline. <i>Adjustment range 95° to 125°. Uses 1 power option.</i>	Standard
---	----------

## STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	FOC			
Active, Coolcore Cover. <i>Solid back, 12 cm contour</i>	Standard			
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	FOC			
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	FOC			
Back Height	50cm. <i>Actual height.</i>	Standard	57cm. <i>Actual height.</i>	FOC

## ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard			
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC			
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	FOC			
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC			
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>	FOC			
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard	
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i>	FOC	
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range 10" to 16". 5"W x 11"L Footplates.</i>	FOC	
Seat pan to heel measurement	(Leave blank for standard factory)	Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard (STD)
--	----------------

## NON-EXPANDABLE CONTROLS

NE + Controls. <i>2 power option ready.</i>	Standard	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

## HAND CONTROL MOUNTS

Fixed.	Standard
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	FOC
Mounting Position.	Left                      FOC                      Right                      Standard

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**MOTABILITY ID 1000791 - 29/09/23 V1**



**Crash Tested**

