



MOTABILITY ID 1000799

Q4 WITH TB FLEX SEAT / POWERED TILT / PRO BACK

DEALERS DETAILS

| | |
|--------------|-----------|
| Dealers Name | |
| Order Date | |
| Address | |
| | Postcode |
| Phone | Email |
| PO No | Order Ref |

DELIVERY DETAILS

| | |
|--------------------|--|
| Dealers Address | |
| Postcode | |
| Quantum Specialist | |
| Clients Weight | |

BASE MODEL

- 7 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors



- ATX Suspension System
- Powered Tilt
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|---------------------------------|----------|
| 3" x 14" Flat Free Black Tyres. | Standard |
| 3" x 14" Pneumatic Black Tyres. | FOC |

COLOURS

| | | | |
|----------------|----------|--|-----|
| Back in Black | Standard | Evolution Orange | FOC |
| Midnight Blue | FOC | Red Corvette | FOC |
| Electric Blue | FOC | Please download our free PQ-365 app to view colours and other product features in more detail. | |
| Deep Purple | FOC | | |
| Pretty in Pink | FOC | | |



SEAT WIDTH

| | | | |
|-----------|-----|-----------|----------|
| 15" Width | FOC | 16" Width | FOC |
| 17" Width | FOC | 18" Width | Standard |
| 19" Width | FOC | 20" Width | FOC |

SEAT DEPTH

| | | | |
|-----------|-----|-----------|----------|
| 15" Depth | FOC | 16" Depth | FOC |
| 17" Depth | FOC | 18" Depth | Standard |
| 19" Depth | FOC | 20" Depth | FOC |

STEALTH SEAT CUSHIONS

| | |
|--|----------|
| No Seat Cushion Required | FOC |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i> | Standard |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

| | |
|---|--|
| Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i> | Standard |
| Seat to Floor Height. <i>From seat pan.</i> | 16" <i>17" Changes Tilt range to -5° to 25°</i> <i>18" Changes Tilt range to -10° to 20°</i> |

STANDARD BACK SUPPORT POSITIONING

| | |
|---|----------|
| Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i> | Standard |
|---|----------|

STEALTH PRO BACK SUPPORTS

| | |
|---|--|
| Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i> | FOC |
| Active, Coolcore Cover. <i>Solid back, 12 cm contour</i> | Standard |
| Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i> | FOC |
| Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i> | FOC |
| Back Height | 50cm. <i>Actual height.</i> Standard 57cm. <i>Actual height.</i> FOC |

ARM SUPPORTS

| | |
|---|---|
| Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i> | Standard |
| TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i> | FOC |
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i> | FOC |
| Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i> | FOC |
| Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i> | FOC |
| Arm Support Pads | 10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC 14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC |

LEG SUPPORTS

| | |
|--|----------|
| 70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i> | Standard |
| Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i> | FOC |
| Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range 10" to 16". 5"W x 11"L Footplates.</i> | FOC |

Seat pan to heel measurement (Leave blank for standard factory) Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

| | |
|--|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | FOC |

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

| | |
|--|----------------|
| Auto Buckle Lap Belt, Push Button Release. | Standard (STD) |
|--|----------------|

NON-EXPANDABLE CONTROLS

| | |
|---|--|
| NE + Controls. <i>2 power option ready.</i> | Standard |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> |
| | Attendant Control. <i>Hand control mounted at rear.</i> |
| | FOC |

HAND CONTROL MOUNTS

| | |
|---|----------|
| Fixed. | Standard |
| Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i> | FOC |
| Mounting Position. | Left |
| | FOC |
| | Right |
| | Standard |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program. <i>Lights and indicators not included.</i> | Standard |
|---|----------|

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 1000799 - 29/09/23 V1



Crash Tested

