

Q6 EDGE 2.0 HD WITH TB3 SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

- 13 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) 4-Pole Motors
- Max User Weight: 28 st 5lbs. (181.0 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|--|----------|
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | Standard |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC |

COLOURS

| | | | |
|---------------|----------|-------------------|-----|
| Back in Black | Standard | Pretty in Pink | FOC |
| Green Machine | FOC | Lemon Crush* | FOC |
| Midnight Blue | FOC | Night Train* | FOC |
| Electric Blue | FOC | Evolution Orange* | FOC |
| Deep Purple | FOC | Military Green* | FOC |
| Red Corvette | FOC | Grey Street* | FOC |

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

20"-24"W

Set width at:

SEAT PACKAGES - SEAT DEPTH

20"-24"D

Set depth at:

STEALTH SEAT CUSHIONS

| | |
|---|-----|
| No Seat Cushion Required | FOC |
| Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i> | STD |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>181.4 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i> | FOC |
| Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>181.4 Kg weight limit.</i> | FOC |
| Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i> | FOC |
| Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>204 Kg weight limit.</i> | FOC |
| Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>226.7 Kg weight limit.</i> | FOC |
| Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

| | |
|--|-----------------------|
| Static Seat. | Standard |
| Static Seat to Floor Height. <i>From seat pan.</i> | 16" 17" 18" |

BACKREST SPECIFICATIONS

| | |
|---|----------|
| TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i> | Standard |
| Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i> | FOC |

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

| | |
|--|----------|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i> | Standard |
|--|----------|

TB3 BACK SUPPORT POSTIONING

| | |
|---|----------|
| TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i> | Standard |
| TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline.</i> | FOC |

BACK CANE ONLY SUPPORT

| | | | | | |
|--|----------|-------|-------|-------|-------|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i> | Standard | | | | |
| Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i> | Standard | | | | |
| Back Cane Height. | 16"H. | 17"H. | 18"H. | 19"H. | 20"H. |
| <i>From seat pan.</i> | 21"H. | 22"H. | 23"H. | 24"H. | 25"H. |

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

| | | | | | |
|--|-------|-----------|-------|-----------|--|
| Backrest Support Width Package. | | 16"-20"W. | | 20"-24"W. | |
| Back Support Width. | 16"W. | 17"W. | 18"W. | 19"W. | |
| | 20"W. | 21"W. | 22"W. | 24"W. | |
| Back Support Height. <i>From seat pan.</i> | 16"H. | 17"H. | 18"H. | 19"H. | |
| | 20"H. | 21"H. | 22"H. | 23"H. | |
| | 24"H. | 25"H. | | | |

19" minimum backrest height when powered recline is selected.

ARM SUPPORTS

| | | | | |
|--|--|----------|--|-----|
| Height Adjustable, Removable, Heavy Duty. | | Standard | | |
| HD Removable Height Package. <i>From seat pan.</i> | 8"-12"H. | FOC | | |
| | 10"-14"H. | Standard | | |
| TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i> | | | FOC | |
| Arm Support Pads | 10" Standard | Standard | 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC |
| | 14" Standard | FOC | 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC |
| | 17" Standard. <i>TB3 only. Minimum seat depth 16".</i> | £50 | | |

CENTRE MOUNTED LEG SUPPORTS

| | | | | |
|---|---|----------|--------------------------------|----------|
| Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with powered lift or tilt.</i> | | Standard | | |
| Manually Adjustable Centre Mount Foot Platform. | | FOC | | |
| Lower Leg Extensions. <i>From seat pan.</i> | 10"-16". | Standard | 13"-19". | FOC |
| Footplates. | 10"W x 8"D Split Footplates. | FOC | 10"W x 11"D Split Footplates. | Standard |
| Calf Support. | Individual Calf Pads. <i>8"W x 6"H.</i> | Standard | Calf Panel. <i>18"W x 8"H.</i> | FOC |

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

| | | | | | | | |
|--|--|---|-----|------------|----------|------------|-----|
| 70° Heavy Duty, Drop-in, Swing-away Leg Supports. <i>Select lowers and footplates below.</i> | | FOC | | | | | |
| 70° HD | Lower Leg Extensions. <i>From seat pan.</i> | 8.5"-12" | FOC | 13"-18" | Standard | 16"-21" | FOC |
| | Footplate Type. | Angle & Depth Adjustable. <i>Select size below.</i> | | | Standard | | |
| | Footplate Size. <i>Angle adj. & Multi-axis only.</i> | 6"W x 8"D | STD | 6"W x 11"D | FOC | 8"W x 11"D | FOC |

INDIVIDUALLY MOUNTED ELEVATING LEG SUPPORTS

| | | | | | | | |
|---|---|-------------|-----|---------|----------|-------------|-----|
| Heavy Duty, Manual Elevating, Swing-away, Leg Supports. <i>Flip-up Extruded footplates 8.5"W x 6"D.</i> | | FOC | | | | | |
| ELEVATING HD | Lower Leg Extensions. <i>From seat pan.</i> | 11"-14" | FOC | 14"-17" | Standard | 12.5"-16.5" | FOC |
| | | 15.5"-19.5" | FOC | | | | |

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

STEALTH HEAD SUPPORTS

| | |
|---|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i> | FOC |
| UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i> | Standard |
| UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ACP150M.</i> | FOC |

Please visit stealthproducts.com for alternative head supports.

PELVIC BELTS

| | |
|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--|----------|

EXPANDABLE CONTROLS

| | |
|--|----------|
| NE + Controls. <i>2 power option ready.</i> | Standard |
| Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |
| Attendant Control. <i>Hand control mounted at rear.</i> | FOC |

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

| | | | | |
|---|----------|-----|-------|----------|
| Fixed. | Standard | | | |
| Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i> | FOC | | | |
| Mounting Position. | Left | FOC | Right | Standard |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program with Lights and indicators. | Standard |
|---|----------|

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
Motability ID 706571 - 10/06/26 V1

