

MOTABILITY ID 706570

## Q6 EDGE 2.0 HD WITH TB3 SEATING / POWER TILT

### DEALERS DETAILS

Dealers Name \_\_\_\_\_  
 Order Date \_\_\_\_\_  
 Address \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Phone \_\_\_\_\_ Email \_\_\_\_\_  
 PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Quantum Specialist \_\_\_\_\_  
 Clients Weight \_\_\_\_\_

### BASE MODEL

- 13 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) 4-Pole Motors
- Max User Weight: 28 st 5lbs. (181.0 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Deep Purple	FOC	Military Green*	FOC
Red Corvette	FOC	Grey Street*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W

Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-24"D

Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>181.4 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>181.4 Kg weight limit.</i>	FOC
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	FOC
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>226.7 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

TRU-Balance 3 HD Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"      18.5"      19.5"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

## TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline.</i>	FOC

## BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
Back Cane Height. <i>From seat pan.</i>	16"H.      17"H.      18"H.      19"H.      20"H. 21"H.      22"H.      23"H.      24"H.      25"H.

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.		16"-20"W.		20"-24"W.	
Back Support Width.	16"W.	17"W.	18"W.	19"W.	
	20"W.	21"W.	22"W.	24"W.	
Back Support Height. <i>From seat pan.</i>	16"H.	17"H.	18"H.	19"H.	
	20"H.	21"H.	22"H.	23"H.	
	24"H.	25"H.			

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

Height Adjustable, Removable, Heavy Duty.		Standard	
HD Removable Height Package. <i>From seat pan.</i>	8"-12"H.	FOC	
	10"-14"H.	Standard	
TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>			FOC
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50	

## CENTRE MOUNTED LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with powered lift or tilt.</i>		Standard	
Manually Adjustable Centre Mount Foot Platform.		FOC	
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	Standard	13"-19". FOC
Footplates.	10"W x 8"D Split Footplates.	FOC	10"W x 11"D Split Footplates. Standard
Calf Support.	Individual Calf Pads. <i>8"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i> FOC

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° Heavy Duty, Drop-in, Swing-away Leg Supports. <i>Select lowers and footplates below.</i>		FOC					
70° HD	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i>				Standard	
	Footplate Size. <i>Angle adj. &amp; Multi-axis only.</i>	6"W x 8"D	STD	6"W x 11"D	FOC	8"W x 11"D	FOC

## INDIVIDUALLY MOUNTED ELEVATING LEG SUPPORTS

Heavy Duty, Manual Elevating, Swing-away, Leg Supports. *Flip-up Extruded footplates 8.5"W x 6"D.* FOC

ELEVATING HD	Lower Leg Extensions.	11"-14"	FOC	14"-17"	Standard	12.5"-16.5"	FOC
	<i>From seat pan.</i>	15.5"-19.5"	FOC				

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. *HMO475-17M / CP10CCFPM.* FOC

UniLink 14" Comfort Plus. *HMO475-17M / CP14CCFPM.* Standard

UniLink 14" Adjustable Comfort Plus. *HMO475-17M / ACP150M.* FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## EXPANDABLE CONTROLS

NE + Controls. *2 power option ready.* Standard

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed. Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.* FOC

Mounting Position.	Left	FOC	Right	Standard
--------------------	------	-----	-------	----------

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators. Standard

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
 Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**Motability ID 706570 - 10/06/24 V1**

