

MOTABILITY ID 714949

ASPEN WITH TB FLEX SEATING / 6 MPH

DEALERS DETAILS

Dealers Name _____
 Order Date _____
 Address _____
 _____ Postcode _____
 Phone _____ Email _____
 PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

 _____ Postcode _____
 Quantum Specialist _____
 Clients Weight _____

BASE MODEL

- 4 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors



- ATX Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 6 Mph Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
--------------------------------	----------

COLOURS

Back in Black	Standard
Electric Blue	FOC
Silver Bullet	FOC
Red Corvette	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified.

SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard
Seat to Floor Height. <i>From seat pan.</i>	17.5" 18" 18.5"

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adjustable range 95° to 125°.</i>	Standard
---	----------

STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
---	----------

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	FOC
Arm Support Pads	10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-away Legrests. <i>Adjustment range 15" to 20".</i>	FOC

Seat pan to heel measurement	(Leave blank for standard factory)	Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.
------------------------------	------------------------------------	--

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.

Standard (STD)

NON-EXPANDABLE CONTROLS

NE + Controls.

Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

HAND CONTROL MOUNTS

Fixed.

Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.*

FOC

Mounting Position. Left FOC Right Standard

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.

Standard

ACCESSORIES

Kerb Climber.

FOC

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITY ID 714949 - 20/07/24 V1

