

R-TRAK WITH TB3 SEATING / POWER TILT

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Front Castors



- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-15"W.
16"-19"W.
20"-22"W.
Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-15"D.
16"-19"D.
20"-22"D.
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	Standard
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.75" 18.75" 19.75"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
--	----------

TB3 BACK SUPPORT POSITIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports, powered recline.</i>	FOC

BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>	Standard
Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	FOC
Back Cane Height. 16"H. 17"H. 20"H.	
<i>From seat pan.</i> 21"H. 24"H. 25"H.	

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.
		18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.	17"H.
	18"H.	19"H.	20"H.	21"H.
	22"H.	23"H.	24"H.	25"H.
	26"H.	27"H.	28"H.	29"H.

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>					Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>					FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>					FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>					FOC
Height Adjustable, Removable, Heavy Duty.					FOC
HD Removable Height Package. <i>From seat pan.</i>		8"-12".	FOC	10"-14".	Standard
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>					FOC
Arm Support Pads.	10" Standard				Standard
	14" Standard				FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>				FOC
	10" Waterfall.				FOC
	14" Waterfall.				FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>					Standard	
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>					FOC	
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". <i>Manually adjustable only.</i>	FOC	13"-19". <i>Manually adjustable only.</i>	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>				FOC	
	10"W x 8"D Split Footplates.				FOC	
	10"W x 11"D Split Footplates.				Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>		FOC	

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>					FOC
90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle & depth Adj. footplates. Adjustment range 7" to 11".</i>					FOC
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>					FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-Away Leg Supports. *Flip-up extruded footplates 8.5"W x 6"D. Select lowers below.* FOC
Minimum seat width 15"

60° TO 90°

Lower Leg Extensions. 11"-13" FOC 13.5"-16.5" Standard 15.5"-19.5" FOC
From seat pan.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. *HMO475-17M / CP08CCFPM.* FOC

UniLink 10" Comfort Plus. *HMO475-17M / CP10CCFPM.* Standard

UniLink 14" Comfort Plus. *HMO475-17M / CP14CCFPM.* FOC

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release. Standard

EXPANDABLE CONTROLS

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. *2 power option ready.* Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.* Standard

Attendant Control. *Hand control mounted at rear.* FOC

HAND CONTROL MOUNTS

Fixed. Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.* FOC

Mounting Position. Left Right

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators. Standard

6 Mph Program with Lights and indicators. FOC


LED Front Castor and Rear Fender Lights. *Includes single switch button mounted to side rail. Not road legal.* Standard

Omit LED Front Castor and Rear Fender Lights. FOC

NOTES & SPECIAL REQUESTS



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
MOTABILITYID 1000800 - 23/07/24 V1



V1: 14/06/2022