

R-TRAK WITH TB3 SEATING / POWER LIFT / POWER TILT

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Front Castors



- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-15"W.

16"-19"W.

20"-22"W.

Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-15"D.

16"-19"D.

20"-22"D.

Set depth at:

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when 10" seat lift is selected.

STEALTH SEAT CUSHIONS

No Seat Cushion Required

FOC

Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. *158.7 Kg weight limit.*

Standard

Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. *160 Kg weight limit.*

FOC

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.*

FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

10" Powered Seat Lift & Tilt. *10" lift & 50° tilt. Uses 2 power options.*

Standard

Powered Lift & Tilt Seat to Floor Height. *From seat pan.*

18.75"

19.75"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. *Select TB3 back options below. Ignore all 'Back Cane' options.*

Standard

Back Canes Only System. *Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.*

FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts.

Standard

Select size in Back Support Dimensions

TB3 BACK SUPPORT POSITIONING

TRU-Balance 3 Static Back. *Angle can be preset manually.*

Standard

TRU-Balance 3 Flip-forward Back. *Not available with TB3 flip-back arm supports, powered recline.*

FOC

BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. *Select size in Back Support Dimensions.*

Standard

Back Cane Only Support Positioning.

Tool Angle Adjustable. *85° to 120° in 5° increments.*

FOC

Back Cane Height.

16"H.

17"H.

20"H.

From seat pan.

21"H.

24"H.

25"H.

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.
		18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.	17"H.
	18"H.	19"H.	20"H.	21"H.
	22"H.	23"H.	24"H.	25"H.
	26"H.	27"H.	28"H.	29"H.

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>					Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>					FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>					FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>					FOC
Height Adjustable, Removable, Heavy Duty.					FOC
HD Removable Height Package. <i>From seat pan.</i>		8"-12".	FOC	10"-14".	Standard
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>					FOC
Arm Support Pads.	10" Standard				Standard
	14" Standard				FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>				FOC
	10" Waterfall.				FOC
	14" Waterfall.				FOC

CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>					Standard	
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>					FOC	
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". <i>Manually adjustable only.</i>	FOC	13"-19". <i>Manually adjustable only.</i>	Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>				FOC	
	10"W x 8"D Split Footplates.				FOC	
	10"W x 11"D Split Footplates.				Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>		FOC	

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>					FOC
90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle & depth Adj. footplates. Adjustment range 7" to 11".</i>					FOC
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>					FOC

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-Away Leg Supports. *Flip-up extruded footplates 8.5"W x 6"D. Select lowers below.*

FOC

Minimum seat width 15"

60° TO 90°

Lower Leg Extensions.
From seat pan.

11"-13"

FOC

13.5"-16.5"

Standard

15.5"-19.5"

FOC

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. *HMO475-17M / CP08CCFPM.*

FOC

UniLink 10" Comfort Plus. *HMO475-17M / CP10CCFPM.*

Standard

UniLink 14" Comfort Plus. *HMO475-17M / CP14CCFPM.*

FOC

PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.

Standard

EXPANDABLE CONTROLS

Q-Logic 3e hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. *2 power option ready.*

Standard

Control Method.

Hand Control. *Select position in Hand Control Mounts below.*

Standard

Attendant Control. *Hand control mounted at rear.*

FOC

HAND CONTROL MOUNTS

Fixed.

Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.*

FOC

Mounting Position.

Left

Right

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.

Standard

6 Mph Program with Lights and indicators.

FOC

LED Front Castor and Rear Fender Lights. *Includes single switch button mounted to side rail. Not road legal.*

Standard

Omit LED Front Castor and Rear Fender Lights.

FOC

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
MOTABILITYID 1000834 - 23/07/24 V1



V1: 14/06/2022