

MOTABILITY ID 707382

## EDGE 3 STRETTO WITH TB3 SEATING / iLEVEL / POWER TILT

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight (required for suspension settings) \_\_\_\_\_

### BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- iLevel Seat Lift and Powered Tilt
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

|  |          |
|--|----------|
| 3" x 12.5" Flat Free Black Tyres.  | Standard |
| 3" x 12.5" Pneumatic Black Tyres.  | FOC      |
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i> | FOC      |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i> | FOC      |

### COLOURS

|               |          |                   |     |
|---------------|----------|-------------------|-----|
| Back in Black | Standard | Pretty in Pink    | FOC |
| Green Machine | FOC      | Lemon Crush*      | FOC |
| Midnight Blue | FOC      | Night Train*      | FOC |
| Electric Blue | FOC      | Evolution Orange* | FOC |
| Blue Jean     | FOC      | Military Green*   | FOC |
| Deep Purple   | FOC      | Grey Street*      | FOC |
| Red Corvette  | FOC      | Raspberry Beret*  | FOC |

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

|               |
|---------------|
| 12"-16"W      |
| 14"-16"W      |
| 16"-20"W      |
| Set width at: |

## SEAT PACKAGES - SEAT DEPTH

|               |
|---------------|
| 12"-16"D      |
| 14"-16"D      |
| 16"-20"D      |
| Set depth at: |

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when iLevel is selected.

## STEALTH SEAT CUSHIONS

|  |          |
|--|----------|
| No Seat Cushion Required   | FOC      |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>                  | Standard |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>                          | FOC      |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC      |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

|  |                       |
|--|-----------------------|
| iLevel Powered Seat Lift & Tilt. <i>12" lift &amp; 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.</i> | Standard              |
| Powered iLevel Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>   | 18"      19"      20" |

## BACKREST SPECIFICATIONS

|   |          |
|---|----------|
| TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>  | Standard |
| Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i> | FOC      |

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a pair of plastic lateral wedge inserts.

## TB3 BACK SUPPORT

|   |          |
|---|----------|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges. <i>Select size in Back Support Dimensions</i> | Standard |
|---|----------|

## TB3 BACK SUPPORT POSITIONING

|  |          |
|--|----------|
| TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>  | Standard |
| TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline. Not suitable for occupied transport.</i> | FOC      |

## BACK CANE ONLY SUPPORT

|   |  |
|---|--|
| Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges. <i>Select size in Back Support Dimensions</i> | Standard                                     |
| Back Cane Only Support Positioning. Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>                               | Standard                                     |
| Back Cane Only Support Positioning. Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>  | FOC  |
| Back Cane Height.   | 16"H      17"H      18"H      19"H      20"H |
| <i>From seat pan.</i>   | 21"H      22"H      23"H      24"H      25"H |

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Ignore Back Support Dimensions if omit back support is selected in TB3 or Back Canes sections. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

| Backrest Support Width Package.               |  | 12"-14"W                                   | 14"-16"W                                   | 16"-20"W                                   |
|---|--|--|--|--|
|   |  | 18"-22"W                                   | 20"-24"W                                   |  |
| Back Support Width.                           | 12"W                                       | 14"W                                       | 16"W                                       | 17"W                                       |
|   | 18"W                                       | 19"W                                       | 20"W                                       | 21"W                                       |
|   | 22"W                                       | 24"W                                       |  |  |
| Back Support Height.<br><i>From seat pan.</i> | 14"H. <i>12"W back only.</i>               | 15"H. <i>12"W back only.</i>               | 16"H.                                      | 17"H.                                      |
|   | 18"H.                                      | 19"H.                                      | 20"H.                                      | 21"H.                                      |
|   | 22"H.                                      | 23"H.                                      | 24"H.                                      | 25"H. <i>Not available with 24"W back.</i> |
|   | 26"H. <i>Not available with 24"W back.</i> | 27"H. <i>Not available with 24"W back.</i> | 28"H. <i>Not available with 24"W back.</i> | 29"H. <i>Not available with 24"W back.</i> |

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

|   |   |          |  |          |
|---|---|----------|--|----------|
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>          |   |          |  | Standard |
| Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>       |   |          |  | FOC      |
| Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i> |   |          |  | FOC      |
| TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>            |   |          |  | FOC      |
| Arm Support Pads  | 10" Standard  | Standard | 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC      |
|   | 14" Standard  | FOC      | 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC      |
|   | 17" Standard. <i>TB3 only. Minimum seat depth 16"</i> | FOC      |  |          |

## CENTRE MOUNTED LEG SUPPORTS

|  |  |          |  |          |
|--|--|----------|--|----------|
| Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i> |  |          |  | Standard |
| Paediatric Adj. Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>   |  |          |  | FOC      |
| Lower Leg Extensions.<br><i>From seat pan.</i>                                 | 8"-12". <i>Paediatric centre mount only.</i>                         | FOC      | 10"-16". <i>Manually adjustable &amp; powered elevating only</i> | Standard |
| Footplates.  | 10"W x 6"D Split Footplates.<br><i>Paediatric centre mount only.</i> | FOC      | 10"W x 8"D Split Footplates.                                     | FOC      |
|  | 10"W x 11"D Split Footplates.  | Standard |  |          |
| Calf Support.  | Individual Calf Pads. <i>6"W x 6"H.</i>                              | Standard |  |          |

## INDIVIDUALLY MOUNTED LEG SUPPORTS

|   |  |  |  |     |
|---|--|--|--|-----|
| 70° Swing-away Legrests. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>                  |  |  |  | FOC |
| 90° Paediatric Swing-away Legrests. <i>4"W x 6"D angle &amp; depth adj. footplates. Adjustment range 7" to 11".</i> |  |  |  | FOC |
| Manual Elevating Swing-Away Legrests. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>        |  |  |  | FOC |

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

|   |         |     |             |          |             |     |
|---|---------|-----|-------------|----------|-------------|-----|
| 60° to 90° Tool Adjustable Swing-away Legrests. <i>Select lowers below. Supplied with 8.5"W x 6"D extruded footplates. Minimum seat width 15"</i> |         |     |             |          |             | FOC |
| 60° to 90° Lower Leg Extensions. <i>From seat pan.</i>  | 11"-13" | FOC | 13.5"-16.5" | Standard | 15.5"-19.5" | FOC |

## STEALTH HEAD SUPPORTS

|  |          |
|--|----------|
| UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i> | Standard |
| UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i> | FOC      |

## PELVIC BELTS

|  |          |
|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--|----------|

## EXPANDABLE CONTROLS

|   |   |  |          |   |   |                      |                     |
|---|---|--|----------|---|---|----------------------|---------------------|
| Q-Logic 3e Controls. <i>3 power option ready.</i>   | Standard  |  |          |   |   |                      |                     |
| <table border="0"> <tr> <td rowspan="2">                     iAccess OLED 4 Button Switch Module.<br/> <i>Mounted to opposite armrest unless specified.<br/>                     Supplied FOC when ilevel is selected.</i> </td> <td rowspan="2">STD</td> <td rowspan="2">}</td> <td>Push Button iAccess</td> <td rowspan="2">If iAccess is selected choose either<br/>Push Button or Toggle lever option.</td> </tr> <tr> <td>Toggle Lever iAccess</td> </tr> </table> | iAccess OLED 4 Button Switch Module.<br><i>Mounted to opposite armrest unless specified.<br/>                     Supplied FOC when ilevel is selected.</i> | STD  | }        | Push Button iAccess                                     | If iAccess is selected choose either<br>Push Button or Toggle lever option. | Toggle Lever iAccess |                     |
| iAccess OLED 4 Button Switch Module.<br><i>Mounted to opposite armrest unless specified.<br/>                     Supplied FOC when ilevel is selected.</i>   |   |  |          | STD   |   | }                    | Push Button iAccess |
|   | Toggle Lever iAccess  |  |          |   |   |                      |                     |
| Omit iAccess OLED 4 Button Switch Module.   |   |  |          |   |   |                      |                     |
| <table border="0"> <tr> <td rowspan="2">Control Method.</td> <td>Hand Control. <i>Select position in Hand Control Mounts below.</i></td> <td>Standard</td> </tr> <tr> <td>Attendant Control. <i>Hand control mounted at rear.</i></td> <td>FOC</td> </tr> </table>  | Control Method.   | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard | Attendant Control. <i>Hand control mounted at rear.</i> | FOC   |                      |                     |
| Control Method.   |   | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |   |   |                      |                     |
|   | Attendant Control. <i>Hand control mounted at rear.</i>   | FOC  |          |   |   |                      |                     |

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

|   |  |       |          |       |          |
|---|--|-------|----------|-------|----------|
| Fixed.  | Standard   |       |          |       |          |
| Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i> | FOC  |       |          |       |          |
| Mounting Position.  | <table border="0"> <tr> <td>Left</td> <td>FOC</td> <td>Right</td> <td>Standard</td> </tr> </table> | Left  | FOC      | Right | Standard |
| Left  | FOC  | Right | Standard |       |          |

## SPEED PROGRAM & LIGHTS PACKAGE

|   |          |
|---|----------|
| 4 Mph Program with Lights and indicators.   | Standard |
| 6 Mph Program with Lights and indicators.   | FOC      |
| LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i> | Standard |
| Omit LED Fender Lights.   | FOC      |

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
 Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**Motability ID 707382 - 21/08/24 V1**

