

Q1450 WITH SYNERGY SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓞ Clients Weight _____

BASE MODEL

- 5 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 3" x 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- 1.5" Anti-tip Wheels
- Swing-arm Suspension System
- VR2 Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Grey Tyres. Standard

REAR CASTORS

2" x 9" Flat Free Grey Tyres. Standard

COLOURS

| | |
|------------------|----------|
| Back in Black | Standard |
| Electric Blue | FOC |
| Red Corvette | FOC |
| Evolution Orange | FOC |
| Pretty in Pink | FOC |

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

| |
|---------------|
| 20"-24"W. |
| 24"-26"W |
| 26"-28"W |
| 28"-30"W |
| 30"-32"W |
| Set width at: |

SEAT PACKAGES - SEAT DEPTH

| |
|---------------|
| 20"-24"D. |
| 24"-28"D |
| Set depth at: |

STEALTH SEAT CUSHIONS

| | |
|---|-----|
| No Seat Cushion Required | FOC |
| Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | STD |
| Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

| | | | | |
|--|-------|-------|-------|----------|
| Static Seat. | | | | Standard |
| Static Seat to Floor Height. <i>From seat pan.</i> | 16.5" | 17.5" | 18.5" | |

SYNERGY BACK SUPPORTS

| | |
|---|----------|
| Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i> | Standard |
| HD Sling Back. <i>Height range 14" to 25" tall from seat pan.</i> | FOC |
| Enter height preference here. | |

BACK SUPPORT POSITIONING

| | | | | | |
|--|----------|-------|-------|-------|-------|
| Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i> | Standard | | | | |
| 1" Endo Cane Mounts. <i>Reduces backrest width by 2" compared to seat width.</i> | FOC | | | | |
| 2" Endo Cane Mounts. <i>Reduces backrest width by 4" compared to seat width.</i> | FOC | | | | |
| Back Cane Height. | 16"H. | 17"H. | 18"H. | 19"H. | 20"H. |
| <i>From seat pan.</i> | 21"H. | 22"H. | 23"H. | 24"H. | 25"H. |

ARM SUPPORTS

| | | | | |
|--|---|----------|----------|----------|
| HD Single Post Tool Height Adjustable Removable. | Standard | | | |
| Bariatric Single Post Tool Height Adjustable Removable. | FOC | | | |
| HD & Bariatric Height Package. <i>From seat pan.</i> | 8"-12". | FOC | 10"-14". | Standard |
| Cantilever Flip-Up Armrests. <i>Adj. range 11"-17". Not available with HD sling back support. Increases overall width by 7".</i> | FOC | | | |
| Arm Support Pads. | 10" Standard. <i>Not available with Bariatric.</i> | FOC | | |
| | 14" Standard. <i>Not available with Bariatric.</i> | Standard | | |
| | 10" Waterfall. <i>Not available with Bariatric.</i> | FOC | | |
| | 14" Waterfall. <i>Not available with Bariatric.</i> | FOC | | |
| | Bariatric. <i>Bariatric arm supports only.</i> | FOC | | |

CENTRE MOUNTED LEG SUPPORTS

17"W x 13"D Base Mounted Flip-up Foot Platform. *Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.* Standard

17"W x 16"D Base Mounted Flip-up Foot Platform. *Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.* FOC

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests. FOC

70° HD Lower Extensions. 8.5"-12". FOC 13"-18". Standard 16"-21". FOC
From seat pan.

70° HD Footplates. 6"W x 8"D. Standard 6"W x 11"D. FOC 8"W x 11"D. FOC

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). *Includes a 60-90° Swing-away leg support hanger for mounting.* LEFT £300 Limb Support £300 RIGHT

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release. Standard

Auto Buckle Lap Belt Length. 60". Standard 70". FOC 80". FOC
90". FOC

NON-EXPANDABLE CONTROLS

VR2 Controls. *With no power options.* Standard

Control Method. Hand Control. *Select position in Hand Control Mounts below.* Standard

HAND CONTROL MOUNTS

Fixed. Standard

Swing-away. *Swings in & out with height, tension and pivot adjustment.* FOC

Mounting Position. Left Right

SPEED PROGRAM

4 Mph Program Standard

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
Motability ID 713665 - 30/07/24 V1



V1: 14/06/2022

