

EDGE 3 STRETTO WITH TB FLEX SEAT / PRO BACK

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight (required for suspension settings) _____

BASE MODEL

- 14 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 12.5" Flat Free Drive Tyres
- 2" x 5" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- NE+ Control System
- 2 x 40 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	FOC
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



STANDARD SEAT DIMENSIONS

15" Width	FOC	16" Width	FOC	15" Depth	FOC	16" Depth	FOC
17" Width	FOC	18" Width	Standard	17" Depth	FOC	18" Depth	Standard
19" Width	FOC	20" Width	FOC	19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard
Manual Seat to Floor Height. <i>From seat pan.</i>	16.5" 17.5" 18.5"

STANDARD BACK SUPPORT POSITIONING *Compatible with standard seat dimensions.*

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°.</i>	Standard
---	----------

STEALTH PRO BACK SUPPORTS *Compatible with standard seat dimensions.*

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	FOC			
Active, Coolcore Cover. <i>Solid back, 12 cm contour</i>	Standard			
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	FOC			
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	FOC			
Back Height	50cm. <i>Actual height.</i>	Standard	57cm. <i>Actual height.</i>	FOC

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard			
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	FOC			
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>	FOC			
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	FOC			
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	FOC			
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
90° Paediatric Swing-away Legrests. <i>Adjustment range 7" to 11".</i>	FOC
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	FOC
Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20".</i>	FOC
Manually Adjustable Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". 5"W x 11"L Footplates.</i>	FOC
Seat pan to heel measurement	(Leave blank for standard factory)

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M/CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M/CP14CCFPM.</i>	FOC

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

NON-EXPANDABLE CONTROLS

NE + Controls. <i>2 power option ready.</i>	Standard
---	----------

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

HAND CONTROL MOUNTS

Fixed.	Standard
--------	----------

Basic swing-away.	FOC
-------------------	-----

Mounting Position.	Left	FOC	Right	Standard
--------------------	------	-----	-------	----------

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
---	----------

NOTES & SPECIAL REQUESTS



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
Motability ID 1000838 - 10/10/24 V1

