

## Q6 EDGE HD WITH SYNERGY SEATING / POWER TILT

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Clients Weight \_\_\_\_\_

### BASE MODEL

- 10 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
35 st 7lbs. (226 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Charcoal Tyres.	Standard
3" x 14" Pneumatic Charcoal Tyres.	FOC

### COLOURS

Back in Black	Standard	Red Corvette	FOC
Night Train	FOC	Pretty in Pink	FOC
Electric Blue	FOC	Evolution Orange	FOC
Midnight Blue	FOC	Grey Street	FOC
Deep Purple	FOC	Green Machine	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W
24"-26"W
26"-28"W
Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-24"D
24"-28"D
Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

HD Powered Seat Tilt. <i>0° to 45°. Uses 1 power option. Only available with R-Net controls.</i>	Standard
HD Tilt Seat to Floor Height. <i>From seat pan.</i>	18.5"      19.5"      20.5"

## SYNERGY BACK SUPPORTS

Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i>	Standard
Enter height preference here.	
Back Cane Height. <i>From seat pan.</i>	14"H.      15"H.      16"H.      17"H.      18"H.      19"H.
	20"H.      21"H.      22"H.      23"H.      24"H.      25"H.

## STANDARD BACK SUPPORT POSTIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
---	----------

## ARM SUPPORTS

HD Single Post Tool Height Adjustable Removable.	Standard
Bariatric Single Post Tool Height Adjustable Removable.	FOC
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12"      FOC      10"-14"      Standard
Arm Support Pads	10" Standard. <i>Not available with Bariatric.</i> FOC      10" Waterfall. <i>Not available with Bariatric.</i> FOC
	14" Standard. <i>Not available with Bariatric.</i> Standard      14" Waterfall. <i>Not available with Bariatric.</i> FOC
	Bariatric. <i>Bariatric arm supports only.</i> FOC

## CENTRE MOUNT LEG SUPPORTS

Adjustable Centre Mount Foot Platform. <i>Mounted to seat.</i>	Standard
Lower Leg Extensions. <i>From seat pan.</i>	10"-16"      Standard      13"-19"      FOC
Footplates.	10"W x 8"D Split Footplates.      FOC      10"W x 11"D Split Footplates.      Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i> FOC      Calf Panel. <i>18"W x 8"H.</i> FOC
	Individual Large Calf Pads. <i>8"W x 8"H tapering to 6"W.</i> Standard

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests.						FOC
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
70° HD Footplates.	6"W x 8"D.	Standard	6"W x 11"D.	FOC	8"W x 11"D.	FOC
HD Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i>						FOC
HD Manual Elevating Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
HD Manual Elevating Footplates.	6"W x 8"D.	FOC	6"W x 11"D.	FOC	8"W x 11"D.	FOC

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard
Auto Buckle Lap Belt Length.	38".	FOC	50".	FOC	60".	Standard
	70".	FOC	80".	FOC	90".	FOC

## EXPANDABLE CONTROLS

R-Net Controls. <i>5 power option ready.</i>		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.		Standard
Basic Swing-away.		FOC
Mounting Position.	Left	Right

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program.	Standard
----------------	----------

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
 Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**Motability ID 721025 - 10/10/24 V1**

