

Technical Bulletin

Bulletin Ref.:

UK-Q2-3-25

Date Issued:

23 May 2025

BULLETIN DETAILS

Product Group: Footplates

Effective Date:

1 May 2025

Model(s) Affected: Powerchairs

Description of Changes/ Diviations:

Quantum Rehab has begun rolling out a new style of footplates across our powerchair range.

Key Improvements

1. Weight Reduction

- The new footplates are 0.6 kg lighter per footplate, reducing overall chair weight.

2. Enhanced Durability

- The previous stick-on coating, which could peel over time, has been replaced with a durable etched coating, more resilient and permanently bonded.



OLD STYLE FOOTPLATES



NEW STYLE FOOTPLATES



New Style Footplate Old Style Footplate

The transition is already in effect. As existing stocks of the old-style footplates run out by size, orders and production are automatically fulfilled with the new version.

Documentation Affected: Owner's Manual Order Form(s) Spec Sheets Other N/A

Link to new version of affected document:

For more detailed information, please contact your local designated Area Sales Manager.

NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

Pride Mobility Products Ltd.
32 Wedgwood Road, Bicester,
Oxfordshire, OX26 4UL

E: sales@quantumrehab.co.uk
T: +44 (0)1869 324600