

## EDGE 3 WITH TB3 SEATING / POWER TILT / POWER RECLINE/ POWER FOOT SUPPORT

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

### BASE MODEL

14 Colour Options  
Height Adj. Removable Arm Supports  
Centre Mount Foot Platform  
TB3 Seating  
14" Flat Free Drive Tyres  
2" x 6" Flat Free Front Castors



2" x 6" Flat Free Rear Castors  
SRS Suspension System  
Q-Logic 3e Control System  
2 x 62 ah Batteries  
4 Mph (6 Km/h) Motors  
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Lemon Crush*	FOC
Midnight Blue	FOC	Night Train*	FOC
Electric Blue	FOC	Evolution Orange*	FOC
Blue Jean	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	Raspberry Beret*	FOC

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

12"-16"W.
14"-16"W.
16"-20"W
20"-22"W
Set width at:

## SEAT PACKAGES - SEAT DEPTH

12"-16"D.
14"-16"D.
16"-20"D
20"-22"D
Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. 0° to 50° COG tilt. Uses 1 power option.				Standard
Powered Tilt Seat to Floor Height. From seat pan.	17.5"	18.5"	19.5"	

## BACK SUPPORT SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below.</i>	Standard
--	----------

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
---	----------

## TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	Standard
---	----------

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Ignore Back Support Dimensions if omit back support is selected in TB3 or Back Canes sections. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.
		18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.		
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.	17"H.
	18"H.	19"H.	20"H.	21"H.
	22"H.	23"H.	24"H.	25"H.
	26"H.	27"H.	28"H.	29"H.
	30"H.	31"H.	32"H.	

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>				FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>				FOC
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				FOC
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with TB3 Flip-back.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with TB3 Flip-back.</i>	FOC

## CENTRE MOUNTED LEG SUPPORTS

Power Elevating Articulating Centre Mount Foot Platform. <i>Uses 1 power option.</i>				FOC
Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	FOC	13"-19".	Standard
Foot Supports.	10"W x 8"D Split Foot Supports.			FOC
	10"W x 11"D Split Foot Supports.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>	FOC

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

Power Elevating Articulating Swing-Away Leg Supports. <i>Includes calf pads.</i>							FOC
ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i>	9"-15"	FOC	13"-21"	Standard		
	Footplate Type.	Angle Adjustable Flip-up. <i>Select size below.</i>					Standard
	Foot Support Size. <i>Includes heel loops.</i>	4"W x 6"D. <i>14" seat widths and above.</i>	FOC	5"W x 7"D. <i>16" seat widths and above.</i>	FOC	6"W x 8"D. <i>18" seat widths and above.</i>	FOC
	Powered Leg Support Operation Method.	Independent. <i>Uses 2 power options.</i>			Combined. <i>Uses 1 power option.</i>		

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	FOC
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	Standard
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	FOC

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC BELTS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

## EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. <i>5 power option ready.</i>		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	FOC

## HAND CONTROL MOUNTS

Fixed.		Standard
Basic Swing-away.		FOC
Mounting Position.	Left	Right

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program with Lights and indicators.	Standard
6 Mph Program with Lights and indicators.	FOC
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

## ACCESSORIES

Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i>	FOC
---	-----

## NOTES & SPECIAL REQUESTS



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**MOTABILITY ID 1003168 - 01/07/25 V1**

