

Technical Bulletin

Bulletin Ref.:

UK-1-4-26

Date Issued:

19 March 2026

BULLETIN DETAILS

Product Group: Quantum Powerchairs

Effective Date:

19 March 2026

Model(s) Affected: All Quantum Powerchairs

Description of Changes/ Diviations:

This bulletin provides information regarding an update to the light only module unit used on the QLNES or Q-logic 3 control system. The unit is being revised to differentiate between a dedicated light only module and a power module with five power options.

Change Details

- The light only module unit will be updated
- The unit will now differ depending on configuration:
 - Light only module
 - 5 power AAM with lights
- This change ensures the correct module is supplied based on the required functionality

Implementation

This change is effective immediately. All new production units will use the updated module configuration, with the correct module supplied based on the specified chair setup.



OLD LIGHT MODULE

PART NUMBER: CTL143986



NEW LIGHT MODULE

PART NUMBER: CTL2005925

Documentation Affected: Owner's Manual Order Form(s) Spec Sheets Other N/A

Link to new version of affected document:

For more detailed information, please contact your local designated Area Sales Manager.

NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

Pride Mobility Products Ltd.
32 Wedgwood Road, Bicester,
Oxfordshire, OX26 4UL

E: sales@quantumrehab.co.uk
T: +44 (0)1869 324600