



# Technical Bulletin

Bulletin Ref.:

UK-2-2-26

Date Issued:

10 April 2026

## BULLETIN DETAILS

Product Group: Applications

Effective Date:

10 April 2026

Model(s) Affected: Econ

### Description of Changes/ Deviations:

With the development of an Android version of Econ, the mobile platform has been rebranded to Econ-M to distinguish it from the Windows-based Econ-W. This update introduces remaining features from Econ-W to Econ-M, along with enhanced security options for easier access.

### KEY CHANGES

- **Platform Rebranding:**
  - Econ-M: New Android mobile app.
  - Econ-W: Remains unchanged for Windows OS.
- **New Features in Econ-M:**
  - Remaining seat programming.
  - Mapped I/O.
  - Monitor: Select and view up to 10 parameters at one time.
  - Battery statistics.
- **System Requirements:**
  - **Apple Devices:** iOS 17 or higher; compatible with iPhone XS/XR, SE 2nd generation, 12 or later; iPad 6th Gen, iPad Air 3rd Gen, iPad Mini 5th Gen, iPad Pro 2nd Gen, or later.
  - **Android Devices:** Android 12 or later.
- **New Features in Econ-W V2:**
  - Aesthetic update to match new Econ mobile
  - Floating window can now be opened with monitor values
  - Floating window can help with setting seating restrictions
- **Security Enhancement:**
  - Ability to use an offline 4-digit PIN code in place of password for login.
  - Upon first download of Econ W V2, log in with username and password, then set the PIN for future use.

Documentation Affected:  Owner's Manual  Order Form(s)  Spec Sheets  Other  N/A

Link to new version of affected document:

*For more detailed information, please contact your local designated Area Sales Manager.*

NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

**Pride Mobility Products Ltd.**  
32 Wedgwood Road, Bicester,  
Oxfordshire, OX26 4UL

E: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
T: +44 (0)1869 324600